

# UNIFIEDGUI

## Release Notes

Benjamin Schiller  
benjamin.bs.schiller@fau.de

April 25, 2022

### **v1.1.0-beta**

**Release date:** 25.04.2022

#### **Minor**

- Decoders: Allow for a restart of the decoder after it has been stopped (all data will be cleared).
- DecoderView: Added a status display for the decoder at the bottom.
- DecoderView: Added a clear button in the decoder toolbar.
- DecoderView/EncoderView: Slight visual changes.
- Decoders: Allow user to specify the x range (default) settings for each decoder individually.
- Plot settings: X range settings are now stored after the program is closed.
- PlotSettingsDialog: Moved step size option to datalines tab.
- PlotView: Slight visual changes.

#### **Bugfixes**

- Fixed a bug when stopping/removing a decoder.

## **v1.0.0-beta**

**Release date:** 21.04.2022

### **Major**

- Encoders: Added encoder functionality — see `Docs/UnifiedGUI.pdf` for detailed information. Added Encoder view to the main view.
- Decoders: Replaced `PARAMETERS` constant by `get.parameters` function.
- Settings: Refactoring of `Settings.py` into `SettingsStore.py` to allow access from all modules without needing a reference to an object.

### **Minor**

- DecoderView: Slight visual changes.
- Live plot: Added the option to adjust the grid of the live plot.
- Live plot: Added the option to adjust the step size of data points drawn in the live plot.
- Menu bar: Added link to release notes PDF to the help menu.

### **Bugfixes**

- DecoderView: Fixed a bug where the subtitle label would not update after removing a decoder.
- DecoderView: Set parameters table to read-only to avoid inconsistencies.

**v1.0.0-alpha**

**Release date:** 03.03.2022